Implant Card.

.....

.....

.....

.....

.....

.....

place the label

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ENGLISH

1. INTRODUCTION

The Implant Card is provided by the Company in paper format to the Health Institutions and it is delivered with the device. The same document is available on the Company website www.limacorporate.com in the menu "USEFUL LINKS" reported in the lower section of the web page.

Any update related to the information provided in the Implant Card and these Instructions is made available to the patient and the Health Institutions through the Company website www.limacorporate.com in the menu "USEFUL LINKS" reported in the lower section of the web page. Updates are also available upon request by e-mail to implantcard@limacorporate.com.

2. INSTRUCTION FOR FILLING

The Implant Card is supplied in a pre-printed format or is downloadable from the Company website; an example is reported in Figure 1.


Figure 1

2.1. PRE-PRINTED PAPER FORMAT

Proceed as follows:

- fill in fields 1, 2, and 3 (see Figure 2);
- apply the specific label delivered with the device (the label indicating "For Implant Card Use Only") in the dedicated space (refer to point 4 of Figure 2);
- fold along the continuous line in the middle of the page layout;
- provide the Implant Card to the patient.

ATTENTION:
If several components are implanted during the surgery, provide an Implant Card for each component, with the proper implant label.

Figure 2

- 1 Patient name or patient ID;
- 2 Name and address of the Health Institution;
- 3 Date of the surgery;
- 4 Device-specific label.

2.2. THE DOWNLOADABLE FORMAT FROM THE COMPANY WEBSITE

Proceed as follows:

- access the following website: www.limacorporate.com;
- access the menu "USEFUL LINKS" reported in the lower section of the web page;
- download and print the Implant Card on paper;
- cut the individual Implant Card along the dotted edges;
- fill in fields 1, 2, and 3 (see Figure 2);
- apply the specific label delivered with the device (the label indicating "For Implant Card Use Only") in the dedicated space (refer to point 4 of Figure 2);
- fold along the continuous line in the middle of the page layout;
- provide the Implant Card to the patient.

ATTENTION: if several components are implanted during the surgery, provide an Implant Card for each component, with the proper implant label.

Health Institutions are required to make available to all patients to whom the device certified according to EC Regulation 2017/745 has been implanted, the following information:

- warnings, precautions, or measures to be taken by the patient or a healthcare professional about reciprocal interference with reasonably foreseeable external influences (e.g. Magnetic resonance imaging - MRI), medical examinations, or environmental conditions. For such information refer to the device Instructions for use;
- any information about the expected lifetime of the device. For such information refer to the device Instructions for use;
- any information about any necessary follow up;
- any other information to ensure the safe use of the device by the patient, including the overall qualitative and quantitative information on materials and substances to which the patient may be exposed. For such information refer to the device Instructions for use.

All the above-mentioned information are available also for devices not certified according to EC Regulation 2017/745, with the exception of the expected lifetime of the device which is provided only for devices certified according to EC Regulation 2017/745.

Health Institutions are responsible for monitoring the Company website for any updates to the above-mentioned additional information.

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The diagram shows a rectangular form titled "Implantaatkaart." with the Lima Corporate logo. It is divided into four numbered sections: 1. A field for patient name or ID, indicated by a person icon and a question mark. 2. A field for the name and address of the care institution, indicated by a person icon with a plus sign. 3. A field for the date of the procedure, indicated by a calendar icon. 4. A large central box labeled "Plaats het etiket." (Place the label). A vertical line runs down the center of the form, and a horizontal line runs across the top, defining the four sections.

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The diagram shows a template for an implant card. It is divided into two main sections by a vertical line. The top left section is titled "Implantaadikaart" and contains three numbered fields: 1 (patient name or ID), 2 (healthcare provider name and address), and 3 (operation date). The top right section features the "LimaCorporate Orthopaedic motion" logo and field 2. The bottom section is a large box labeled "pange kleebis siia" (stick the sticker here) with field 4 (device location sticker) at the bottom center. A vertical dashed line runs down the middle, and a horizontal dashed line runs across the top of the bottom section.

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The diagram shows a rectangular form titled "Implantatkort" on the left. On the right side, there is the Lima Corporate logo with the tagline "Orthopaedic motion". The form is divided into four numbered sections:

- 1:** A field for patient name or ID, indicated by a person icon and a question mark.
- 2:** A field for the care provider's name and address, indicated by a person icon with a plus sign.
- 3:** A field for the date of the operation, indicated by a calendar icon.
- 4:** A large rectangular area at the bottom labeled "placera etiketten" (place the label).

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The diagram shows a form titled "Karta implantu" with the Lima Corporate logo. It is divided into four numbered sections: 1. Patient name/ID, 2. Facility name and address, 3. Date of procedure, and 4. A large box for the device label. A vertical dashed line runs down the center, and a horizontal dashed line runs across the top. A large number 4 is placed in a circle at the bottom center of the label box.

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The diagram shows a rectangular implant card template with a dashed border. At the top left, it says "implanto kortelė". At the top right is the "Lima Corporate Orthopaedic motion" logo. The card is divided into four numbered sections: 1. A field for patient name or number, indicated by a person icon and a question mark. 2. A field for the orthopedic institution name and address, indicated by a person icon with a plus sign. 3. A field for the operation date, indicated by a calendar icon. 4. A large rectangular area at the bottom for pasting a label, with the instruction "užklijuokite etiketę" (paste the label). A vertical dashed line runs down the center of the card, and a horizontal dashed line runs across the middle, intersecting at the top of section 4. The number "2 pav." is written at the bottom right of the diagram.

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The diagram shows the front of an Implant Card with the following fields:

- Field 1:** A person icon followed by a line for the patient's name or ID.
- Field 2:** A person icon with a plus sign followed by a line for the name of the institution.
- Field 3:** A calendar icon followed by a line for the date of surgery.
- Field 4:** A large rectangular area labeled "poġġi t-tikketta" (place the sticker) at the bottom, intended for a specific device sticker.

The card also features the title "Tessera tal-Impjant" and the "Lima Corporate Orthopaedic motion" logo.

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The diagram shows a rectangular form titled "Implanttikortti" (Implant Card) with the Lima Corporate logo. It is divided into four numbered sections:

- 1:** A field for patient name or ID, indicated by a person icon and a question mark.
- 2:** A field for the name and address of the healthcare provider, indicated by a person icon with a plus sign.
- 3:** A field for the surgery date, indicated by a calendar icon.
- 4:** A large rectangular area at the bottom labeled "Liimaa tarra" (Stick the label), intended for a "For Implant Card Use Only" label.

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Limacorporate S.p.A.